

Model Name	Model #	Material	Design	Injector	Cartridge	A-Constant Optical	Diopter Range	HCPCS Code
<b>Presbyopia Correcting</b>								
Tecnis Symphony OptiBlue	DXR00V	Hydrophobic Acrylic	1-piece	Preloaded		119.3	5.0 - 34.0 *0.5 D steps	V2788 Presbyopia correcting function of IOL
Tecnis Symphony OptiBlue Toric II	DXW 150,225,300,375	Hydrophobic Acrylic	1-piece	Preloaded		119.3	5.0 - 34.0 *0.5 D steps	V2788 Presbyopia correcting function of IOL
Tecnis Synergy	DFR00V	Hydrophobic Acrylic	1-piece	Preloaded		119.3	5.0 - 34.0 *0.5 D steps	V2788 Presbyopia correcting function of IOL
Tecnis Synergy Toric II	DFW 150,225,300,375	Hydrophobic Acrylic	1-piece	Preloaded		119.3	5.0 - 34.0 *0.5 D steps	V2788 Presbyopia correcting function of IOL
Tecnis Multifocal 2.75 add	ZKB00	Hydrophobic Acrylic	1-piece	Platinum DK7796	1MTEC30	119.3	5.0 - 34.0 *0.5 D steps	V2788 Presbyopia correcting function of IOL
Tecnis Multifocal Toric II 2.75 add	ZKU 150,225,300,375	Hydrophobic Acrylic	1-piece	Platinum DK7796	1MTEC30	119.3	5.0 - 34.0 *0.5 D steps	V2788 Presbyopia correcting function of IOL
Tecnis Multifocal 3.25 add	ZLB00	Hydrophobic Acrylic	1-piece	Platinum DK7796	1MTEC30	119.3	5.0 - 34.0 *0.5 D steps	V2788 Presbyopia correcting function of IOL
Tecnis Multifocal Toric II 3.25 add	ZLU 150,225,300,375	Hydrophobic Acrylic	1-piece	Platinum DK7796	1MTEC30	119.3	5.0 - 34.0 *0.5 D steps	V2788 Presbyopia correcting function of IOL
<b>Toric - Monofocal</b>								
Tecnis Eyhance Toric II	DIU 150,225,300,375, 450,525,600	Hydrophobic Acrylic	1-piece	Preloaded		119.3	5.0 - 34.0 *0.5 D steps	V2787 Astigmatism correcting function of an IOL□
Tecnis Toric II	ZCU 150,225,300,375, 450,525,600	Hydrophobic Acrylic	1-piece	Platinum DK7796	1MTEC30	119.3	5.0 - 34.0 *0.5 D steps	V2787 Astigmatism correcting function of an IOL□
<b>Monofocal</b>								
Tecnis Eyhance	DIB00	Hydrophobic Acrylic	1-piece	Preloaded		119.3	5.0 - 34.0 *0.5 D steps	V2632 Posterior chamber IOL
Tecnis Simplicity	DCB00	Hydrophobic Acrylic	1-piece	Preloaded		119.3	5.0 - 34.0 .5 D steps	V2632 Posterior chamber IOL
Tecnis One	ZCB00	Hydrophobic Acrylic	1-piece	Platinum DK7796	1MTEC30	119.3	5.0 - 34.0 *0.5 D steps	V2632 Posterior chamber IOL
Tecnis Acrylic	ZA9003	Hydrophobic Acrylic	3-piece	EMERALDXL	Emerald C (Green Box)	119.1	10.0 - 30.0 *0.5 D steps	V2632 Posterior chamber IOL
Sensar1	AAB00	Hydrophobic Acrylic	1-piece	DK7796	1MTEC30	119.0	6.0 - 30.0 *0.5 D steps	V2632 Posterior chamber IOL
Sensar Sensar (low diopter) Sensar (low diopter)	AR40e AR40E AR40M	Hydrophobic Acrylic	3-piece	EMERALDXL	Emerald C (Green Box)	118.7	6.0 - 30.0 *0.5 D steps +2.0 - +5.5 *0.5 D steps -10.0 - +1.5 * 0.5 D steps	V2632 Posterior chamber IOL